



Freeze fracture setup

<https://search.labfacilities.wur.nl/SearchDetail.aspx?deviceid=9ecbd1ae-782a-4d90-b335-723a896dab5d>

Brand

Leica

Type

EM BAF060

Contact

Marcel Giesbers (marcel.giesbers@wur.nl)

Organisation

Plant Sciences

Department

Wageningen Electron Microscopy Centre

Description

The Leica EM BAF060 is a high-end sample preparation system for freeze fracturing, etching, and freeze drying of double replica, high resolution carbon/metal mix coatings for TEM and SEM analysis.